

Protect and preserve cells and

tissues at low temperatures



Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands

T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com



Product info

ROKEPIE[®] = a bioregulator to preserve cells & tissue at 2-8° C.

ROKEPIE® Admirated of 9726 GNGM Art.nr: RS.01.10 Product: ROKEPIE® Charge: 16H24 Exp.date: 01-dec-17 RESEARCH ONLY STORE AT 2-8 °C FOR M

Additive

 \circ No need to change medium

- Non-toxic
 - \odot No need to remove before rewarming (37°C
 - Supports cell growth after rewarming
- RUO product
 - For Research Use Only
- Fully synthetic
 - No animal components, antibiotics and preservatives
 - One defined, small molecule: SUL-109
- Sterile additive
 - Open, Add, Incubate & Use
- Non-GMP



Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com



Application



Create stock solution





Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com



Use and benefits





Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com



Use and benefits





Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands

T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com



Hypothermia rewarming model

- Assess cell protective effect
- Read-outs:
 - Trypan blue exclusion (dead cells)
 - Neutral red uptake (viable cells)
 - Growth curve after storage at 2-8 °C





Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands

T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com



Tested primary cell lines

HUVEC control cells stored in 96-wells plate for 3 days at 2-8° C.







HUES9 control cells after storage in 6-wells plate for 7 days at 2-8° C.



HUES9 cells after storage in 6-wells plate for 7 days at 2-8° C with ROKEPIE[®].





Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands

T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com



Tested cell lines (2)

hADSC control cells stored for 11 days at 2-8° C.



hADSCs stored for 11 days at 2-8° C with ROKEPIE[®].



NRK 52E control cells after storage for 21 days at 2-8° C.



NRK 52E cells after storage for 21 days at 2-8° C with ROKEPIE[®].



ROKEPIE[®]

Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com



Overview test results

	Concentration	Storage duration (days)
Primary astrocytes	1:10	3
CaCo-2	1:10/1:100	7
HEK 293	1:10/1:100	7
HUVEC	1:10/1:100	13
MS5	1:10	7
N2A*	1:10/1:100	7
	0.1/1/10 μΜ	
NRK 52E	1:10/1:100	21
SH-Sy5y	1:10/1:50	7
hADSC	1:100	3
HDF*	1:20	2
THP-1*	1:25	2
Human skin punches*	1:10	5
Aorta rings rat*	100 μΜ	N.A.
Pigs heart*	10 mL/L	1

ROKEPIE[®]

Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com



Interesting topics

- Cell culture media

 Other than DMEM, RPMI, Ham's F-12
- Culture media composition

 Serum, buffers
- Suspension cell culture
 THP-1
- Functionality of cell lines

 Differentiation, cytokine production
- Detachment of intact monolayer
 As seen with HUVEC
- Passage number and cold tolerance



Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands

T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com



Background: hibernation technology

Natural resistance against damage

- Hibernation is nature's solution for surviving harsh times (food, temperature)
- Hibernators use a protective mechanism during hibernation

Interesting phenomena of hibernation

- No impact of fattening or cooling on health
- No signs of damage to organs or brain
- Active inhibition of metabolism strong decrease in food and oxygen consumption
- Various adaptations to protect body and mind



ROKEPIE[®]

Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands

T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com

Users about ROKEPIE

European Society of Gene & Cell Therapy

"ROKEPIE <u>effectively promotes</u> cell viability of primary human hematopoietic stem cells during hypothermic storage"





Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands

T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com





Optimization of methods for hypothermic preservation of primary human Hematopoietic Stem Cells in culture

Aysun Sarikaya¹, Trui Visser¹, Özgür Özyüncü², Duygu Uçkan-Çetinkaya^{1,3}, Gerard Wagemaker¹, Fatima Aerts-Kaya¹

Results

- Center for Stem Cell Research and Development, Department of Stem Cell Sciences, Graduate School of Health Sciences, Hacettepe University, Ankara, Turkey Department of Obstetrics and Gynecology, Faculty of Medicine, Hacettepe University, Ankara, Turkey
- Pediatric Bone Marrow Transplantation Unit, Division of Pediatric Hematology, Ihsan Dogramaci Children's Hospital, Faculty of Medicine, Hacettepe University, Ankara, Turkey 3)
- nding author e-mail: fatimaaerts@yahoo.com

Introduction

Biogenic amines (dopamine and serotonin) protect cells against hypothermia-induced damage (1). The newly developed compound ROKEPIE®-S01 (Sulfateq BV) has been demonstrated to be successful in hypothermic storage and transportation of several cell lines, as well as thrombocytes (2), but not of primary cells, including hematopoietic stem cells (HSCs). Hypothermic preservation of HSCs could be preferred over short-term cryopreservation, to prevent cell loss during freezing/thawing and could be particularly useful for shortterm storage, e.g., during conditioning of transplant patients or transport of HSC transplants. Here, we assessed optimal doses of this compound and culture conditions to preserve CD34+ umbilical cord blood (UCB) HSCs.



Figure 1. Prolonged hypothermic storage of Hematopoletis Stem Cells To optimize conditions for prolonged hypothermic storage of HSCs, umbilical cord blood CD34+ cells were selected and cultured in presence of different growth factor cocktalis, cultured for 2 or 7 days, supplemented without washing steps with Rokeple at dilutions of 1:5 – 1:50 and analyzed for colony forming units, CD34+ content, apoptosis (Annexin-V-FITC/Propidium lodide) and cell cycle status.

Materials and Methods

CD34+ HSCs were isolated from freshly collected umbilical cord blood samples using magnetically labeled anti-CD34 antibodies (Miltenyi). HSCs were supplemented with TPO, TPO/SCF, TPO/SCF/Flt3L or TPO/SCF/aFGF/IGF2/Angptl3/heparin (STIFA3) for 2 days (analogous to our gene transfer method) or 7 days (analogous to our HSC expansion method). After 2 or 7 day cell culture, cell counts, colony forming assays (Methocult, Stem Cell Technologies) and CD34 cell numbers were assessed using a FACSARIA (Becton Dickinson). ROKEPIE®-S01 was then added directly to the medium, the plates closed airtight and stored for 4-7 days at 4°C. After rewarming, the frequency of apoptotic cells and cell cycle status were measured (Figure 1)

Acknowledgements

This study was supported by grants from the EU 7th Framework program CELL-PID no 261387, The Turkish Research Council TÜBİTAK 2221 no 2015/8 and the Turkish Governement DPT K-120 640: PEDI-STEM

ROKEPIE[®]

Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands

T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com



ROKEPIE® is registered by Sulfateg BV

Two day cultures with TPO resulted in the lowest frequency of apoptotic cells after hypothermic storage in the presence of 1:10 or 1:20 diluted ROKEPIE@S01. After 7 days of culture/hypothermic storage, cell viability was 10-30% higher in cultures containing the compound, whereas numbers of cycling cells were 10-20% lower, in particular when stimulated with STIFA3 (Figure 2).



Figure 2. Apoptosis in HSC cultures after hypothermic storage CD34+ UCB-CD34+ cells were cultured in presence of TPO for 2 days or STIFA3 for 7 days. ROKEPIE was added to the cultures and the flasks were closed air-tight and stored at 4°C for 3 or 7 days. Apoptosis was measured using Annexin-V-FITC and Propidium lodide

In short-term HSC cultures, followed by short-term storage at 4°C, ROKEPIE@S01 does not appear to affect cell viability or cell cycle state. Under long-term HSC expansion conditions, followed by 7 days of hypothermic storage, however, effects of the compound are significant, as evident from decreased rate of apoptosis, and increased levels of cells in G0/G1 phase (Table 1).

Table 1.	Cell c	ycle analy	sis of HSC	cultures after	hypothermic storage

Cell cycle	99/91 (%)	\$ (%)	02/M (%)	\$H92M (%)
2 days TPO/3 days 4°C without Rokepie	12,2	11,7	11,8	23,5
2 days TPO/3 days 4°C Rokepie	15,1	15,2	15,0	90,2
2 day TPO/7 days 4°C without Rokapia	19,9	11,5	8,5	20,0
2 day TPO/7 days 4 [®] C Rokepie	14,5	15,9	16,4	32,3
7 days STIFA3/3 days 4°C without Rokepie	30,5	27,3	23,9	51,2
7 days STIFA3/3 days 4*C Rokepie	27,7	28,5	26,6	53,1
7 days STIFA3/7 days 4*C without Rokeple	12,7	25,1	31,4	58,5
7 days STIFA37 days 4*C Rokepie	24,6	26,1	24,3	50,4

Conclusions

ROKEPIE®-S01 effectively promotes cell viability of primary HSCs during hypothermic storage, and is currently assessed for its effects on long-term repopulating HSC.

Rerences

Talaei F et al., PLoS One. 2011;6(7):e22568 Talaei F et al. Patent application WO2011128458 A1

Users about ROKEPIE

"We like to share our encouraging results; When ROKEPIE[®] is used at 1/10 dilution and dispensed in both apical and basal wells of the plate, we see a major improvement. Without adding ROKEPIE[®] we simply cannot store and ship in cells in our gel formulation."

"We tested the ROKEPIE[®] with primary human keratinocytes, for 3 days at 4°C. We found that the cells recommenced normal proliferation after they were returned to the incubator. It is of interest for transport purposes, which we will investigate further."

"Thus far, I am <u>quite satisfied</u> with the results that we have been achieving using the ROKEPIE[®] hibernation reagent (especially also given the ease of its use) and I'd be happy to perform a few more experiments. "

"We have tested ROKEPIE[®] in a 1:20 dilution added to our cell culture by storing them in the fridge for 24 hrs. The morphology after storage was <u>very good</u>, many mitoses were visible. The 3D model images after 10 and 14 days were very good, as well as the barrier function."



Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com



Want more information about ROKEPIE? Please visit our website or contact us through info@rokepie.com

Hibernation applied in life sciences

WWW.ROKEPIE.COM



Admiraal de Ruyterlaan 5 9725 GN Groningen The Netherlands

T: +31 50 313 7464 E: info@rokepie.com W: www.rokepie.com

